

# Viscotester VT-06

# Uni-Cylinder Rotational Viscotester Allows Easy Measurement of Fluid Viscosity



The VT-06 is designed for quality control applications in the manufacturing process of industrial products such as petrochemicals, paint, and adhesives, as well as foodstuffs. Viscosity measurements covering a wide range are possible, such as gear oil used in construction machinery. Measurement is performed by simply submerging a rotor in the fluid. The resistance to rotor movement caused by the viscosity (torque) is measured to obtain direct readings.

- Compact and light weight make the unit easily portable and allow operation with one hand
- Can be powered by alkaline batteries, nickel-hydride rechargeable batteries, or AC adapter
- Direct indication of viscosity in decipascal-seconds (SI units)
- Dedicated stand for measurement available as option

# **(Usage)**

- 1. Attach rotor to unit and hold unit in the hand or place on dedicated stand. (Unit should be approximately horizontal in either case.)
- 2. Insert rotor in sample fluid, turn power on, and select rotor number.
- 3. Press start button and read indicated viscosity.
- \*The supplied extension rod can reach fluid that is further away. (Only for use with the No.1 and No.2 rotors.)

### Specifications

•			
Measurement range	No. 3 rotor: 0.3 to 13 dPa·s (with No. 3 cup)		
	No. 1 rotor: 3 to 150 dPa·s (with JIS 300 mL beaker*1)		
	No. 2 rotor: 100 to 4000 dPa·s (with JIS 300 mL beaker*1)		
Sample fluid capacity	No. 1 and No. 2 rotor	(with JIS 300 mL beaker*1)	approx. 300 mL
	No. 3 rotor	(with No. 3 cup)	approx. 150 mL
	Clearance between rotor end and cup bottom:		
		about 15 mm	
Measurement accuracy	±10 %±1 digit of indicated value, reproducibility ±5 %		
Rotor speed	62.5 rpm		
Power supply	IEC LR6 (size AA) alkaline batteries,		
	nickel-hydride rechargeable batteries, AC adapter VA-05JA		

Dimensions and Weight	175 (H) × 77 (W) × 40 (D) mm (without protruding parts),		
	Approx. 260 g (without batteries)		
Supplied accessories	No. 1 rotor (dia. 24 × 53 × 166 mm)	SUS304	1
	No. 2 rotor (dia. 15 x 1 x 113 mm)	SUS304	1
	No. 3 rotor (dia. 45 × 47 × 160 mm)	SUS304	1
	No. 3 Cup (dia. 52.6 x 75 mm)	SUS304	1
	Extension rod (900 mm · 300× 3)	SUS304	1
	IEC LR6 (size AA) alkaline batteries		4


Note: The Viscotester cannot measure accurately with anything other than supplied cups or the JIS 300 mL beaker.

\*1 JIS R 3503 : 1994, φ78×103 (H)

#### Options

Product name	Product number  VA-04  VA-05JA	
Stand		
AC adapter		

## Rotors and Cups (unit: mm)



#### Sample amount for measurement

No.3 Cup	approx.150 mL		
Commercially available 300 mL beaker	approx.350 mL		

Note: For certain fluids, readings may differ slightly from other viscometers, depending on properties of target fluids, mechanical factors, as well as specific gravity, rotor speed, and other aspects.

# ■ Viscotester measurement examples (for reference)

Product type	Viscosity	Rotor		
Newtonian fluids				
Castor oil	6 dPa·s	No.3		
Starch syrup	1000 dPa·s	No.2		
Non-Newtonian fluids				
Condensed milk	16 dPa·s	No.1		
Chocolate syrup	25 dPa·s	No.1		
Tomato ketchup	43 dPa·s	No.1		
Pure honey	76 dPa·s	No.1		
Toothpaste	320 dPa·s	No.2		
Starch paste	310 dPa·s	No.2		

 $<sup>*\,\</sup>mbox{Measurement}$  temperature: 23 °C

# ■CGS Unit and SI Unit

$$1cP = \frac{1}{1.000}Pa \cdot s = 0.01 dPa \cdot s$$
  $1P = \frac{1}{10}Pa \cdot s = 1 dPa \cdot s$ 

P(poise), cP(centi poise), Pa·s(pascal-seconds), dPa·s(decipascal-seconds)



JCSS JCSS 0197 RION CO., LTD. is recognized by the JCSS which uses ISO/IEC 17025 as an accreditation standard and bases its accreditation scheme on ISO/IEC 17011. JCSS is operated by the accreditation body (IA Japan) which is a signatory to the Asia Pacific Accreditation Cooperation (APAC) as well as the International Laboratory Accreditation Cooperation (ILAC). The Quality Assurance Section of RION CO., LTD. is an international MRA compliant JCSS operator with the accreditation number JCSS 0197.



\* Specifications subject to change without notice





3-20-41, Higashimotomachi, Kokubunji, Tokyo 185-8533, Japan Tel: +81-42-359-7888 Fax: +81-42-359-7442