

**LABORHITS**

**Mini-3D Rocker PS-3D**

**Space-saving**

**Robust**

**Quiet**



**Features**

- Simple speed adjustment
- Robust drive
- Suitable for continuous operation up to 7 days
- Low energy consumption
- With anti-slip mat
- Fixed angle of 7°

Platform shaker with smooth, 3-dimensional tumbling movement. Ideal for mixing tubes, culture bottles, culture dishes and microtiter plates.

**Specifications**

**Type of movement:** tumbling  
**Tumbling angle:** 7°  
**Speed range:** 5 ... 60 rpm  
**Max. Load:** 1.0 kg  
**Ambient Temperature:** 4 ... 40 °C  
**Dimensions Plattform (W x D):** 215 x 215 mm  
**Dimensions (W x D x H):** 85 x 150 x 80 mm  
**Weight:** 1.2 kg  
**Power supply:** 12 V

| Type                   | PK | Cat. No.  |
|------------------------|----|-----------|
| PS-3D, with EU/UK-plug | 1  | 9.721 014 |