

LABORHITS

VITgrip™ lab bottles, PP, GL 45

Features

- Break resistant
- Slim, waisted shape facilitates the handling
- Molded volume graduations support slip-safety while working with gloves
- Leakproof** due to optimized sealing system of the bottle thread and the screw cap with ergonomic design
- Easy readability of the volume at all times due to molded, double-sided graduations with high accuracy ($\pm 5\%$)
- Smooth pouring of liquids and easy to clean due to hydrophobic material and rounded shape without dirt catching edges
- With tamper-evident closure



Excellent handling

Leakproof

** VITgrip™ lab bottles are leakproof (tested with distilled water at room temperature and atmospheric pressure and in combination with the supplied VITLAB® screw closure).



Allround lab bottle made of polypropylene with GL 45 thread.

Very good chemical resistance against most acids, bases and alcoholic solutions. Furthermore, the VITgrip™ lab bottles and the supplied VITLAB® tamper-evident closures (PP) are suitable for contact with foodstuffs (according to regulation EU No. 10/2011).

Capacity ml	Height* mm	Bottom diam. mm	PK	Cat. No.
125	103	54	6	6.270 627
250	149	64	6	6.270 628
500	192	77	6	6.270 629
1000	234	97	6	6.270 630
2000	278	126	1	6.270 631

* Height without screw cap