

LABORHITS

Sampling Dippers DispoDipper LaboPlast® / SteriPlast®

**Sampling and transport
in one device**

Safe handling



NEW

Features

- PP, blue transparent
- Individually packaged and optionally sterilised by gamma rays (SteriPlast®)
- Production and packaging according to clean room class 7
- Complies with EU food requirements and FDA regulations
- Scaled sampling tube up to 50 ml
- Immersion depth 270 mm

Sampling tube manufactured in a clean room specially designed for disposable use. Due to the blue colour, it is also perfect for use in the food, animal feed and pharmaceutical industries. After sampling, the container can be closed and is leakproof. The sampling tube is suitable for 2nd drum openings and therefore can be used for sampling barrels, canisters and other containers.



DispoDipper is completely manufactured in one piece and can be used for sampling and subsequent transport of the sample in one device. In this way, cross-contamination or sample impurities are almost impossible. After sampling, the handle at the predetermined breaking point can be detached from the sample container. This allows for easy and space-saving transport of the sample.

Description	Capacity ml	Diam. mm	Height mm	Overall length mm	PK	Cat. No.
LaboPlast®	50	40	98	417	20	6.290 524
SteriPlast®	50	40	98	417	20	6.290 525