

# LABORHITS

## Duran® YOUTILITY Bottle System

**Complete system with purging ring  
and GL 45 screw cap**

**Specially shaped gripping  
zones for safer handling**

**Made of DURAN®  
borosilicate 3.3 glass**



# LABORHITS

## Duran® YOUTILITY Bottle System


Specially shaped gripping zones on both sides of the bottle for easier and safer handling. With the DURAN® YOUTILITY bottle thread opening or closing the bottle is significantly faster. The thread is fully compatible with DIN GL 45 closures and other accessories. The slimmer DURAN® YOUTILITY bottle shape allows more optimal use of limited space in autoclaves and laboratory refrigerators. A pre-defined labelling area is compatible with the dedicated DURAN® self-adhesive DURAN® YOUTILITY bottle labels. Nominal volume is shown at the top of the easy-to-read graduation scale for fast determination of the volumes. Glass bottle body is moulded from the tried and tested DURAN® borosilicate 3.3 Type 1 neutral glass.

DURAN® glass offers very good chemical resistance and high temperature resistance.

The DURAN® YOUTILITY bottle has an integrated Retrace Code to allow full product traceability. The lot specific certificate of conformance is available at: [www.duran-group.com](http://www.duran-group.com).

Bottles are available in either clear or amber coloured glass. The amber colour meets the requirements of the USP <660> and Eur. Ph. requirements for UV light protection.

Each DURAN® YOUTILITY bottle is supplied as a complete system, with a pouring ring (cyan polypropylene) and a GL 45 cap (cyan polypropylene).



Capacity ml	Colour	Ext. diam mm	Height mm	PK	Cat. No.
125	clear	55	124	4	<b>9.071 679</b>
250	clear	66	158	4	<b>9.071 678</b>
500	clear	78	193	4	<b>9.071 677</b>
1000	clear	93	253	4	<b>9.071 676</b>
125	amber	55	124	4	<b>9.071 672</b>
250	amber	66	158	4	<b>9.071 673</b>
500	amber	78	193	4	<b>9.071 674</b>
1000	amber	93	253	4	<b>9.071 675</b>

### Accessories

Description	PK	Cat. NO.
Spare screw caps, GL 45, cyan, PP, food safe	10	<b>9.209 137</b>
Spare pouring ring, GL 45, cyan, PP	16	<b>9.209 138</b>
Bottle tags, GL 45, 8 colours (2 tags each), silicone	16	<b>9.209 139</b>
Self-adhesive label set, 200 bottle cap labels, PES	1	<b>9.209 140</b>

