

LABORHITS

Overhead Mixer LLG-uniROTATOR 2



**Very Flexible
Clamp Fixation System**

Easy to Clean

Cost-Efficient

Also order

LLG-Centrifuge Tubes *ECONOMY*
(6.270 403, 6.270 404, 6.270 405)

Analogue overhead mixer with vertical rotary motion for gentle but effective mixing. Ideal for biological samples. Suitable for cutting and mixing processes of powdered and liquid samples in tubes up to 50 ml. With easily interchangeable clamps that allow a quick cleaning of the mixer.

- Compact design
- Housing easy to clean
- Non-slip, adjustable screw feet

Scope of delivery:

Overhead mixer incl. 12 clamps each
for 2 ml, 15 ml and 50 ml tubes

Specifications

Type of movement: overhead

Variable speed: 0 - 80 rpm

Rotation angle: 360°

Permissible shaking weight (incl. attachment): 2 kg

Motor rating input: 20 W

Dimensions (W x D x H): 530 x 150 x 190 mm

Power supply: 100 - 240 V, 50/60 Hz

Protection class: IP21

Type	PK	Cat. No.
LLG-uniROTATOR 2 with EU-plug	1	6.263 660
LLG-uniROTATOR 2 with UK-plug	1	6.263 661
Clamp for 2 ml centrifuge tube	12	6.263 662
Clamp for 15 ml centrifuge tube	12	6.263 663
Clamp for 50 ml centrifuge tube	12	6.263 664