

## 1 Acidimeter model Jozzi

Steam current distiller, complete with glassware and reagents, for volatile and total acidity detection. According to EEC norms.

<u>Technical features</u>: Heating system: electric Power: 400 W

According to EEC regulations

Description	PK	Cat. No.
Jozzi 1-place	1	A113421



## 2 Accessori e ricambi per acidimetri Jozzi

Description	PK	Cat. No.
Glass boiler	1	X200731
Glass boiling tube	1	X200732
Flask 250 ml, with spherical joint	1	X200736
Beaker 250 ml, for distillation collection	1	X200737
Cooling element, spherical joints	1	X200075
Junction for multi-position Jozzi, having electrical heating system	1	X200076
Junction for 1-place Jozzi, having electrical heating system	1	X200097
Clamp for spherical joint 29/15	1	X200070
Clamp for spherical joint 18/9	1	X200088
Heating element for 1-place Jozzi	1	X200232
Heating element for multi-position Jozzi	1	X200738

Boiler, tube, cooler, junction

